



SPLIT TUBE SYSTEM

18 mm & 22 mm
DIAMETERS




 **New Pivoting Frame Up
To 40 Degrees Of Distal
Blade Rotation**

BASIC SPLIT TUBE SET

Basic Set Includes


18 & 22mm Diameters,
50, 60, 70, 80mm Lengths



EXTENDED TUBE SET

Additional Items:

90, 100, 110, 120mm Lengths.
Stepdown Frame, Locking
Block, 3rd Arm Attachment,
Hooks & Blades.



Optional Items

Articulating Arm



Light Cable Ports Offer Superior Directional Light

